

SONICLEAN™

SO EASY - SO FAST - SO CLEAN




超声波数码桌上型清洗机

超声波数码桌上型清洗机适用于净洁医疗和牙医手术工具。备有1.6公升至80公升容量大小，以及密封结构和可编程功能供选择。

Ultrasonic Digital Bench Top Machines

Ultrasonic digital benchtop machines are suitable for cleaning medical and dental surgical instruments. Comes in sizes from 1.6 litre to 80 litre capacities.



喷淋式超声波清洗机

喷淋式超声波清洗机采用超声波，配合将干净清洁液灌入工具中冲洗的组合方法，用于清洗中空套管插入工具以及腹腔镜检工具。“内外同时净洁”

Ultrasonic Irrigator Machine

Ultrasonic irrigator used for cleaning cannulated hollow instruments and laparoscopic instruments. "Cleans the inside as well as the outside"



脉冲扫描功率

ADAPTIVE PULSE SWEEP POWER

SONICLEAN™
SO EASY - SO FAST - SO CLEAN



超声波起重清洗机

超声波起重清洗机是用于清理沉重的手术工具，包括用于进行锁孔手术的腹腔镜检工具。它也具备双槽装置供冲洗和烘干。

Ultrasonic Lifter Washer Machine

Ultrasonic Lifter Washer machine used for cleaning heavy loads of surgical instruments including laparoscopic instruments used for key hole surgery. Also comes in double tank configuration for rinsing and drying.

超声波性能测量仪

超声波性能测量仪用于测量仪器设备的超声波活动以验证其效率，比其他测试方法更快速廉宜。



Ultrasonic Performance Meter


Ultrasonic Performance Meter (UPM) used to measure the performance of ultrasonic activity in the machine to validate efficiency. Faster and cheaper than alternative testing methods.

欲知有关上述产品详情，请上网浏览：
www.soniclean.com.au

或请联络：

Soniclean Pty Ltd
ABN 89 057 730 917
38 Anderson Street
Thebarton SA 5031
Australia

电话: +61 8 8234 8398
传真: +61 8 8234 8391
sales@soniclean.com.au
www.soniclean.com.au



TGA Registration
Class 1 Certification
ARTG # 104498

PRODUCT OF AUSTRALIA